

Antibodies Product Datasheet

Product Name:	DAND5 Antibody
Catalog Number:	C50100
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human DAND5
Concentration:	1 mg/ml
Gene Name:	DAND5 CER2 CKTSF1B3 GREM3 SP1
Protein Name:	DAN domain family member 5 (Cerberus-like protein 2) (Cerl-2) (Cysteine knot superfamily 1, BMP antagonist 3) (Gremlin-3)
Other Name:	
Human Swiss Prot No.:	Q8N907
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of DAND5 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year