

Antibodies Product Datasheet

Product Name:	HMGCLL1 Antibody
Catalog Number:	C50141
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human HMGCLL1
Concentration:	1 mg/ml
Gene Name:	HMGCLL1
Protein Name:	3-hydroxymethyl-3-methylglutaryl-CoA lyase, cytoplasmic (EC 4.1.3.4) (3-hydroxymethyl-3-methylglutaryl-CoA lyase-like protein 1) (Endoplasmic reticulum 3-hydroxymethyl-3-methylglutaryl-CoA lyase) (er-cHL)
Other Name:	
Human Swiss Prot No.:	Q8TB92
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of HMGCLL1 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year