

Antibodies Product Datasheet

Product Name:	BHLHA15 Antibody
Catalog Number:	C50163
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human BHLHA15
Concentration:	1 mg/ml
Gene Name:	BHLHA15 BHLHB8 MIST1
Protein Name:	Class A basic helix-loop-helix protein 15 (bHLHa15) (Class B basic helix-loop-helix protein 8) (bHLHb8) (Muscle, intestine and stomach expression 1) (MIST-1)
Other Name:	
Human Swiss Prot No.:	Q7RTS1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of BHLHA15 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year