

Antibodies Product Datasheet

Product Name:	HYAL4 Antibody
Catalog Number:	C50176
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human HYAL4
Concentration:	1 mg/ml
Gene Name:	HYAL4
Protein Name:	Hyaluronidase-4 (Hyal-4) (EC 3.2.1.35) (Chondroitin sulfate endo-beta-N-acetylgalactosaminidase) (Chondroitin sulfate hydrolase) (CSHY) (Hyaluronoglucosaminidase-4)
Other Name:	
Human Swiss Prot No.:	Q2M3T9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of HYAL4 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year