

Antibodies Product Datasheet

Product Name:	ASGR2 Antibody
Catalog Number:	C50180
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ASGR2
Concentration:	1 mg/ml
Gene Name:	ASGR2 CLEC4H2
Protein Name:	Asialoglycoprotein receptor 2 (ASGP-R 2) (ASGPR 2) (C-type lectin domain family 4 member H2) (Hepatic lectin H2) (HL-2)
Other Name:	
Human Swiss Prot No.:	P07307
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ASGR2 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year