

Antibodies Product Datasheet

Product Name:	KIF14 Antibody
Catalog Number:	C50181
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human KIF14
Concentration:	1 mg/ml
Gene Name:	KIF14 KIAA0042
Protein Name:	Kinesin-like protein KIF14
Other Name:	
Human Swiss Prot No.:	Q15058
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of KIF14 at Human
Storage and Stability:	Store at -20° C/1 year