

Antibodies Product Datasheet

Product Name:	ALG6 Antibody
Catalog Number:	C50208
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ALG6
Concentration:	1 mg/ml
Gene Name:	ALG6 My046
Protein Name:	Dolichyl pyrophosphate Man9GlcNAc2 alpha-1,3-glucosyltransferase (EC 2.4.1.267) (Asparagine-linked glycosylation protein 6 homolog) (Dol-P-Glc:Man(9)GlcNAc(2)-PP-Dol alpha-1,3-glucosyltransferase) (Dolichyl-P-Glc:Man9GlcNAc2-PP-dolichyl glucosyltransferas
Other Name:	
Human Swiss Prot No.:	Q9Y672
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ALG6 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year