

Antibodies Product Datasheet

Product Name:	BRF2 Antibody
Catalog Number:	C50214
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human BRF2
Concentration:	1 mg/ml
Gene Name:	BRF2 BRFU PRO1470
Protein Name:	Transcription factor IIIB 50 kDa subunit (TFIIIB50) (hTFIIIB50) (B-related factor 2) (BRF-2) (hBRFU)
Other Name:	
Human Swiss Prot No.:	Q9HAW0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of BRF2 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year