

Antibodies Product Datasheet

Product Name:	GLYAT Antibody
Catalog Number:	C50252
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human GLYAT
Concentration:	1 mg/ml
Gene Name:	GLYAT ACGNAT CAT GAT
Protein Name:	Glycine N-acyltransferase (EC 2.3.1.13) (EC 2.3.1.71) (Acyl-CoA:glycine N-acyltransferase) (AAG) (Aralkyl acyl-CoA N-acyltransferase) (Aralkyl acyl-CoA:amino acid N-acyltransferase) (Benzoyl-coenzyme A:glycine N-acyltransferase) (Glycine N-benzoyltransferase)
Other Name:	
Human Swiss Prot No.:	Q6IB77
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of GLYAT at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year