

Antibodies Product Datasheet

Product Name:	IER3IP1 Antibody
Catalog Number:	C50258
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human IER3IP1
Concentration:	1 mg/ml
Gene Name:	IER3IP1 HSPC039
Protein Name:	Immediate early response 3-interacting protein 1
Other Name:	
Human Swiss Prot No.:	Q9Y5U9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of IER3IP1 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year