

Antibodies Product Datasheet

Product Name:	KIRREL3 Antibody
Catalog Number:	C50261
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human KIRREL3
Concentration:	1 mg/ml
Gene Name:	KIRREL3 KIAA1867 NEPH2 UNQ5923/PRO4502/PRO19814
Protein Name:	Kin of IRRE-like protein 3 (Kin of irregular chiasm-like protein 3) (Nephrin-like protein 2)
Other Name:	
Human Swiss Prot No.:	Q8IZU9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of KIRREL3 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year