

Antibodies Product Datasheet

Product Name:	LPPR3 Antibody
Catalog Number:	C50264
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human LPPR3
Concentration:	1 mg/ml
Gene Name:	LPPR3 PHP2 PRG2
Protein Name:	Lipid phosphate phosphatase-related protein type 3 (EC 3.1.3.4) (PAP-2-like protein 2) (Plasticity-related gene 2 protein) (PRG-2)
Other Name:	
Human Swiss Prot No.:	Q6T4P5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of LPPR3 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year