

Antibodies Product Datasheet

Product Name:	Cyclon/CCDC86 Antibody
Catalog Number:	C50356
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human Cyclon/CCDC86
Concentration:	1 mg/ml
Gene Name:	CCDC86 CYCLON
Protein Name:	Coiled-coil domain-containing protein 86 (Cytokine-induced protein with coiled-coil domain)
Other Name:	
Human Swiss Prot No.:	Q9H6F5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Cyclon/CCDC86 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year