

Antibodies Product Datasheet

Product Name:	IKIP Antibody
Catalog Number:	C50383
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human IKIP
Concentration:	1 mg/ml
Gene Name:	IKBIP IKIP
Protein Name:	Inhibitor of nuclear factor kappa-B kinase-interacting protein (I kappa-B kinase-interacting protein) (IKBKB-interacting protein) (IKK-interacting protein)
Other Name:	
Human Swiss Prot No.:	Q70UQ0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of IKIP at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year