

Antibodies Product Datasheet

Product Name:	ART5 Antibody
Catalog Number:	C50391
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ART5
Concentration:	1 mg/ml
Gene Name:	ART5 UNQ575/PRO1137
Protein Name:	Ecto-ADP-ribosyltransferase 5 (EC 2.4.2.31) (ADP-ribosyltransferase C2 and C3 toxin-like 5) (ARTC5) (Mono(ADP-ribosyl)transferase 5) (NAD(P)(+)-arginine ADP-ribosyltransferase 5)
Other Name:	
Human Swiss Prot No.:	Q96L15
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ART5 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year