

Antibodies Product Datasheet

Product Name:	DNAJC4 Antibody
Catalog Number:	C50401
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human DNAJC4
Concentration:	1 mg/ml
Gene Name:	DNAJC4 HSPF2 MCG18
Protein Name:	DnaJ homolog subfamily C member 4 (DnaJ-like protein HSPF2) (Multiple endocrine neoplasia type 1 candidate protein number 18)
Other Name:	
Human Swiss Prot No.:	Q9NNZ3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of DNAJC4 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year