

Antibodies Product Datasheet

Product Name:	FICD Antibody
Catalog Number:	C50404
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human FICD
Concentration:	1 mg/ml
Gene Name:	FICD HIP13 HYPE UNQ3041/PRO9857
Protein Name:	Adenosine monophosphate-protein transferase FICD (EC 2.7.7.n1) (AMPylator FICD) (FIC domain-containing protein) (Huntingtin yeast partner E) (Huntingtin-interacting protein 13) (HIP-13) (Huntingtin-interacting protein E)
Other Name:	
Human Swiss Prot No.:	Q9BVA6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of FICD at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year