

Antibodies Product Datasheet

Product Name:	KDSR Antibody
Catalog Number:	C50411
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human KDSR
Concentration:	1 mg/ml
Gene Name:	KDSR FVT1
Protein Name:	3-ketodihydrosphingosine reductase (KDS reductase) (EC 1.1.1.102) (3-dehydrosphinganine reductase) (Follicular variant translocation protein 1) (FVT-1)
Other Name:	
Human Swiss Prot No.:	Q06136
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of KDSR at Human, Mouse
Storage and Stability:	Store at -20° C/1 year