

Antibodies Product Datasheet

Product Name:	DPP8 Antibody
Catalog Number:	C50448
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human DPP8
Concentration:	1 mg/ml
Gene Name:	DPP8 DPRP1 MSTP097 MSTP135 MSTP141
Protein Name:	Dipeptidyl peptidase 8 (DP8) (EC 3.4.14.5) (Dipeptidyl peptidase IV-related protein 1) (DPRP-1) (Dipeptidyl peptidase VIII) (DPP VIII) (Prolyl dipeptidase DPP8)
Other Name:	
Human Swiss Prot No.:	Q6V1X1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of DPP8 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year