

Antibodies Product Datasheet

Product Name:	KTELC1 Antibody
Catalog Number:	C50452
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human KTELC1
Concentration:	1 mg/ml
Gene Name:	POGLUT1 C3orf9 CLP46 KTELC1 MDSRP MDS010 UNQ490/PRO1006
Protein Name:	Protein O-glucosyltransferase 1 (EC 2.4.1.-) (CAP10-like 46 kDa protein) (hCLP46) (KTEL motif-containing protein 1) (Myelodysplastic syndromes relative protein)
Other Name:	
Human Swiss Prot No.:	Q8NBL1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of KTELC1 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year