

Antibodies Product Datasheet

Product Name:	C18orf21 Antibody
Catalog Number:	C50461
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human C18orf21
Concentration:	1 mg/ml
Gene Name:	C18orf21 XTP13 PNAS-124 PNAS-131
Protein Name:	UPF0711 protein C18orf21 (HBV X-transactivated gene 13 protein) (HBV XAg-transactivated protein 13)
Other Name:	
Human Swiss Prot No.:	Q32NC0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of C18orf21 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year