

Antibodies Product Datasheet

Product Name:	C19orf36 Antibody
Catalog Number:	C50462
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human C19orf36
Concentration:	1 mg/ml
Gene Name:	IZUMO4 C19orf36 UNQ831/PRO1758
Protein Name:	Izumo sperm-egg fusion protein 4 (Sperm 22 kDa protein c113)
Other Name:	
Human Swiss Prot No.:	Q1ZYL8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of C19orf36 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year