

Antibodies Product Datasheet

Product Name:	KHDC1 Antibody
Catalog Number:	C50468
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human KHDC1
Concentration:	1 mg/ml
Gene Name:	KHDC1 C6orf148
Protein Name:	KH homology domain-containing protein 1
Other Name:	
Human Swiss Prot No.:	Q4VXA5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of KHDC1 at Human
Storage and Stability:	Store at -20° C/1 year