

Antibodies Product Datasheet

Product Name:	ANKLE1 Antibody
Catalog Number:	C50472
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ANKLE1
Concentration:	1 mg/ml
Gene Name:	ANKLE1 ANKRD41 LEM3
Protein Name:	Ankyrin repeat and LEM domain-containing protein 1 (Ankyrin repeat domain-containing protein 41) (LEM-domain containing protein 3)
Other Name:	
Human Swiss Prot No.:	Q8NAG6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ANKLE1 at Human
Storage and Stability:	Store at -20° C/1 year