

Antibodies Product Datasheet

Product Name:	ABCE1 Antibody
Catalog Number:	C50480
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ABCE1
Concentration:	1 mg/ml
Gene Name:	ABCE1 RLI RNASEL1 RNASELI RNS4I OK/SW-cl.40
Protein Name:	ATP-binding cassette sub-family E member 1 (2'-5'-oligoadenylate-binding protein) (HuHP68) (RNase L inhibitor) (Ribonuclease 4 inhibitor) (RNS4I)
Other Name:	
Human Swiss Prot No.:	P61221
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ABCE1 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year