

Antibodies Product Datasheet

Product Name:	ARL6IP6 Antibody
Catalog Number:	C50481
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ARL6IP6
Concentration:	1 mg/ml
Gene Name:	ARL6IP6 PFAAP1
Protein Name:	ADP-ribosylation factor-like protein 6-interacting protein 6 (ARL-6-interacting protein 6) (Aip-6) (Phosphonoformate immuno-associated protein 1)
Other Name:	
Human Swiss Prot No.:	Q8N6S5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ARL6IP6 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year