

Antibodies Product Datasheet

Product Name:	FBXO44 Antibody
Catalog Number:	C50493
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human FBXO44
Concentration:	1 mg/ml
Gene Name:	FBXO44 FBG3 FBX30 FBX44 FBX6A FBXO6A
Protein Name:	F-box only protein 44 (F-box protein FBX30) (F-box/G-domain protein 3)
Other Name:	
Human Swiss Prot No.:	Q9H4M3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of FBXO44 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year