

Antibodies Product Datasheet

Product Name:	GNGT2 Antibody
Catalog Number:	C50495
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human GNGT2
Concentration:	1 mg/ml
Gene Name:	GNGT2 GNG8 GNG9 GNGT8
Protein Name:	Guanine nucleotide-binding protein G(I)/G(S)/G(O) subunit gamma-T2 (G gamma-C) (G-gamma-8) (G-gamma-9) (Guanine nucleotide binding protein gamma transducing activity polypeptide 2)
Other Name:	
Human Swiss Prot No.:	O14610
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of GNGT2 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year