

Antibodies Product Datasheet

Product Name:	ISCA2 Antibody
Catalog Number:	C50528
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ISCA2
Concentration:	1 mg/ml
Gene Name:	ISCA2 HBLD1
Protein Name:	Iron-sulfur cluster assembly 2 homolog, mitochondrial (HESB-like domain-containing protein 1)
Other Name:	
Human Swiss Prot No.:	Q86U28
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ISCA2 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year