

Antibodies Product Datasheet

Product Name:	JP45/JSRP1 Antibody
Catalog Number:	C50529
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human JP45/JSRP1
Concentration:	1 mg/ml
Gene Name:	JSRP1 JP45
Protein Name:	Junctional sarcoplasmic reticulum protein 1 (Junctional-face membrane protein of 45 kDa homolog) (JP-45)
Other Name:	
Human Swiss Prot No.:	Q96MG2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of JP45/JSRP1 at Human
Storage and Stability:	Store at -20° C/1 year