

## Antibodies Product Datasheet

|                               |                                                                                                                            |
|-------------------------------|----------------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | GAB3 Antibody                                                                                                              |
| <b>Catalog Number:</b>        | C50534                                                                                                                     |
| <b>Reactivity:</b>            | Human, Mouse                                                                                                               |
| <b>Applications:</b>          | WB                                                                                                                         |
| <b>Dilutions:</b>             | WB 1:500-2000                                                                                                              |
| <b>Modification:</b>          | Unmodified/Total                                                                                                           |
| <b>Source:</b>                | Polyclonal, Rabbit,IgG                                                                                                     |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.      |
| <b>Immunogen:</b>             | Synthesized peptide derived from human GAB3                                                                                |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                                    |
| <b>Gene Name:</b>             | GAB3                                                                                                                       |
| <b>Protein Name:</b>          | GRB2-associated-binding protein 3 (GRB2-associated binder 3) (Growth factor receptor bound protein 2-associated protein 3) |
| <b>Other Name:</b>            |                                                                                                                            |
| <b>Human Swiss Prot No.:</b>  | Q8WWW8                                                                                                                     |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                    |
| <b>Specificity:</b>           | This antibody detects endogenous levels of GAB3 at Human, Mouse                                                            |
| <b>Storage and Stability:</b> | Store at -20° C/1 year                                                                                                     |