

Antibodies Product Datasheet

Product Name:	ARAP2 Antibody
Catalog Number:	C50535
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ARAP2
Concentration:	1 mg/ml
Gene Name:	ARAP2 CENTD1 KIAA0580
Protein Name:	Arf-GAP with Rho-GAP domain, ANK repeat and PH domain-containing protein 2 (Centaurin-delta-1) (Cnt-d1) (Protein PARX)
Other Name:	
Human Swiss Prot No.:	Q8WZ64
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ARAP2 at Human
Storage and Stability:	Store at -20° C/1 year