

Antibodies Product Datasheet

Product Name:	HEYL Antibody
Catalog Number:	C50543
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human HEYL
Concentration:	1 mg/ml
Gene Name:	HEYL BHLHB33 HRT3
Protein Name:	Hairy/enhancer-of-split related with YRPW motif-like protein (hHeyL) (Class B basic helix-loop-helix protein 33) (bHLHb33) (Hairy-related transcription factor 3) (HRT-3) (hHRT3)
Other Name:	
Human Swiss Prot No.:	Q9NQ87
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of HEYL at Human, Mouse
Storage and Stability:	Store at -20° C/1 year