

Antibodies Product Datasheet

Product Name:	CBR4 Antibody
Catalog Number:	C50547
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human CBR4
Concentration:	1 mg/ml
Gene Name:	CBR4
Protein Name:	Carbonyl reductase family member 4 (EC 1.-.-) (3-oxoacyl-[acyl-carrier-protein] reductase) (EC 1.1.1.-) (Quinone reductase CBR4)
Other Name:	
Human Swiss Prot No.:	Q8N4T8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CBR4 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year