

Antibodies Product Datasheet

Product Name:	FTSJ2 Antibody
Catalog Number:	C50555
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human FTSJ2
Concentration:	1 mg/ml
Gene Name:	FTSJ2 FJH1
Protein Name:	Putative ribosomal RNA methyltransferase 2 (EC 2.1.1.-) (Protein ftsJ homolog 2) (rRNA (uridine-2'-O-)-methyltransferase)
Other Name:	
Human Swiss Prot No.:	Q9UI43
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of FTSJ2 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year