

Antibodies Product Datasheet

Product Name:	ABI3BP Antibody
Catalog Number:	C50574
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ABI3BP
Concentration:	1 mg/ml
Gene Name:	ABI3BP NESHBP TARSH
Protein Name:	Target of Nesh-SH3 (Tarsh) (ABI gene family member 3-binding protein) (Nesh-binding protein) (NeshBP)
Other Name:	
Human Swiss Prot No.:	Q7Z7G0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ABI3BP at Human
Storage and Stability:	Store at -20° C/1 year