

Antibodies Product Datasheet

Product Name:	DAP3 Antibody
Catalog Number:	C50592
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human DAP3
Concentration:	1 mg/ml
Gene Name:	DAP3 MRPS29
Protein Name:	28S ribosomal protein S29, mitochondrial (MRP-S29) (S29mt) (Death-associated protein 3) (DAP-3) (Ionizing radiation resistance conferring protein)
Other Name:	
Human Swiss Prot No.:	P51398
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of DAP3 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year