

Antibodies Product Datasheet

Product Name:	ENTPD8 Antibody
Catalog Number:	C50597
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ENTPD8
Concentration:	1 mg/ml
Gene Name:	ENTPD8 UNQ2492/PRO5779
Protein Name:	Ectonucleoside triphosphate diphosphohydrolase 8 (E-NTPDase 8) (NTPDase 8) (NTPDase8) (EC 3.6.1.5)
Other Name:	
Human Swiss Prot No.:	Q5MY95
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ENTPD8 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year