

Antibodies Product Datasheet

Product Name:	ICMT Antibody
Catalog Number:	C50605
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ICMT
Concentration:	1 mg/ml
Gene Name:	ICMT PCCMT
Protein Name:	Protein-S-isoprenylcysteine O-methyltransferase (EC 2.1.1.100) (Isoprenylcysteine carboxymethyltransferase) (Prenylated protein carboxyl methyltransferase) (PPMT) (Prenylcysteine carboxyl methyltransferase) (pcCMT)
Other Name:	
Human Swiss Prot No.:	O60725
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ICMT at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year