

Antibodies Product Datasheet

Product Name:	LARGE Antibody
Catalog Number:	C50611
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human LARGE
Concentration:	1 mg/ml
Gene Name:	LARGE KIAA0609 LARGE1
Protein Name:	Glycosyltransferase-like protein LARGE1 (EC 2.4.-.-) (Acetylglucosaminyltransferase-like 1A)
Other Name:	
Human Swiss Prot No.:	O95461
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of LARGE at Human, Mouse
Storage and Stability:	Store at -20° C/1 year