

Antibodies Product Datasheet

Product Name:	METTTL11A Antibody
Catalog Number:	C50618
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human METTTL11A
Concentration:	1 mg/ml
Gene Name:	NTMT1 C9orf32 METTTL11A NRMT AD-003
Protein Name:	N-terminal Xaa-Pro-Lys N-methyltransferase 1 (EC 2.1.1.244) (Alpha N-terminal protein methyltransferase 1A) (Methyltransferase-like protein 11A) (N-terminal RCC1 methyltransferase) (X-Pro-Lys N-terminal protein methyltransferase 1A) (NTM1A)
Other Name:	
Human Swiss Prot No.:	Q9BV86
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of METTTL11A at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year