

Antibodies Product Datasheet

Product Name:	HMGN4 Antibody
Catalog Number:	C50695
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human HMGN4
Concentration:	1 mg/ml
Gene Name:	HMGN4 HMG17L3 NHC
Protein Name:	High mobility group nucleosome-binding domain-containing protein 4 (Non-histone chromosomal protein HMG-17-like 3) (Non-histone chromosomal protein)
Other Name:	
Human Swiss Prot No.:	O00479
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of HMGN4 at Human
Storage and Stability:	Store at -20° C/1 year