

Antibodies Product Datasheet

Product Name:	EFTUD1 Antibody
Catalog Number:	C50719
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human EFTUD1
Concentration:	1 mg/ml
Gene Name:	EFTUD1 EFL1 FAM42A
Protein Name:	Elongation factor Tu GTP-binding domain-containing protein 1 (Elongation factor-like 1) (Protein FAM42A)
Other Name:	
Human Swiss Prot No.:	Q7Z2Z2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of EFTUD1 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year