

Antibodies Product Datasheet

Product Name:	GAGE2D Antibody
Catalog Number:	C50756
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human GAGE2D
Concentration:	1 mg/ml
Gene Name:	GAGE2D; GAGE8
Protein Name:	G antigen 2D (GAGE-2D) (Cancer/testis antigen 4.8) (CT4.8) (G antigen 8) (GAGE-8)
Other Name:	
Human Swiss Prot No.:	Q9UEU5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of GAGE2D at Human
Storage and Stability:	Store at -20° C/1 year