

Antibodies Product Datasheet

Product Name:	BORIS Antibody
Catalog Number:	C50763
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human BORIS
Concentration:	1 mg/ml
Gene Name:	CTCFL BORIS
Protein Name:	Transcriptional repressor CTCFL (Brother of the regulator of imprinted sites) (CCCTC-binding factor) (CTCF paralog) (CTCF-like protein) (Cancer/testis antigen 27) (CT27) (Zinc finger protein CTCF-T)
Other Name:	
Human Swiss Prot No.:	Q8NI51
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of BORIS at Human, Mouse
Storage and Stability:	Store at -20° C/1 year