

Antibodies Product Datasheet

Product Name:	KIAA0774 Antibody
Catalog Number:	C50767
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human KIAA0774
Concentration:	1 mg/ml
Gene Name:	MTUS2 CAZIP KIAA0774 TIP150
Protein Name:	Microtubule-associated tumor suppressor candidate 2 (Cardiac zipper protein) (Microtubule plus-end tracking protein TIP150) (Tracking protein of 150 kDa)
Other Name:	
Human Swiss Prot No.:	Q5JR59
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of KIAA0774 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year