

Antibodies Product Datasheet

Product Name:	LCE1B Antibody
Catalog Number:	C50769
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human LCE1B
Concentration:	1 mg/ml
Gene Name:	LCE1B LEP2 SPRL2A
Protein Name:	Late cornified envelope protein 1B (Late envelope protein 2) (Small proline-rich-like epidermal differentiation complex protein 2A)
Other Name:	
Human Swiss Prot No.:	Q5T7P3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of LCE1B at Human
Storage and Stability:	Store at -20° C/1 year