

Antibodies Product Datasheet

Product Name:	ATP6V0A4 Antibody
Catalog Number:	C50780
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ATP6V0A4
Concentration:	1 mg/ml
Gene Name:	ATP6V0A4 ATP6N1B ATP6N2
Protein Name:	V-type proton ATPase 116 kDa subunit a isoform 4 (V-ATPase 116 kDa isoform a4) (Vacuolar proton translocating ATPase 116 kDa subunit a isoform 4) (Vacuolar proton translocating ATPase 116 kDa subunit a kidney isoform)
Other Name:	
Human Swiss Prot No.:	Q9HBG4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ATP6V0A4 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year