

Antibodies Product Datasheet

Product Name:	ANXA9 Antibody
Catalog Number:	C50785
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ANXA9
Concentration:	1 mg/ml
Gene Name:	ANXA9 ANX31
Protein Name:	Annexin A9 (Annexin XXXI) (Annexin-31) (Annexin-9) (Pemphaxin)
Other Name:	
Human Swiss Prot No.:	O76027
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ANXA9 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year