

Antibodies Product Datasheet

Product Name:	DMRT1 Antibody
Catalog Number:	C50788
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human DMRT1. AA range:330-373
Concentration:	1 mg/ml
Gene Name:	DMRT1 DMT1
Protein Name:	Doublesex- and mab-3-related transcription factor 1 (DM domain expressed in testis protein 1)
Other Name:	
Human Swiss Prot No.:	Q9Y5R6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of DMRT1 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year